

BERCHTOLD

Podo



PODIATRY DRILL WITH
BRUSHLESS VACUUM TECHNOLOGY



Berchtold Podo Q Podiatry Drill
with Dust Extraction

- 25,000 rpm Handpiece
- Handpiece On/Off Switch
- LCD Touch Screen Display
- Clockwise/Counter-clockwise Rotation
- Brushless Vacuum Motor Technology
- Touch Button Controls
- 3 Programmable Memory Buttons for Suction Level, Burr Speed & Burr Direction
- Handpiece Holder

- Height 8cm
- Width 33cm (includes h/p holder)
- Depth 25cm
- Weight 2.4kg
- Handpiece Weight 140g

Suction Performance:

190 L/min

Noise Level (full vacuum):

70-75 dB

Warranty:

3 years

Manufactured in:

Germany

